



RESIN 8 COLORS

A SUPERIOR RESIN INSTRUMENT


How the best perform



**YOUR PERFORMANCE
IS OUR PASSION.**

At Hu-Friedy, every instrument, from its proprietary metal and textured design, to the color for organization, is designed with one purpose—to help you perform at your best.

When we created Resin 8 Colors scalers and curettes, we combined resin handles in 8 eye-catching colors with features that maximize your comfort and efficiency to achieve the best clinical outcomes—every patient, every day.



RESIN 8 COLORS POINTS OF PERFORMANCE



EVEREDGE® TECHNOLOGY

Stays sharper longer, which means less time sharpening

STEEL INNER CORE DESIGN

For added strength and tactile sensitivity

LARGE DIAMETER RESIN HANDLE

Optimum weight for a more comfortable grip and reduced pinch force

UNIQUE TEXTURED HANDLE DESIGN

Wave grooves and unique knurling texture increase rotational control and provide a light secure grasp

SIGNATURE SERIES® GRIPS

Silicone grip provides more comfort and color identification

PERFECT ANGLES

Consistently hand crafted blades for the most efficient calculus removal

COLOR CODED SYSTEM

To easily identify scalers and curettes that correspond to anatomical areas of the mouth



RESIN 8 COLORS GRACEY COLOR CODING

Resin 8 Colors Gracey curettes are organized by areas of the mouth, providing easy identification.



GRACEYS: **1/2, 3/4**

1/2 and 3/4 Graceys available in Standard, Rigid, After Five®, After Five Rigid, Mini Five® and Mini Five Rigid varieties



GRACEYS: **5/6**

5/6 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid varieties



GRACEYS: **7/8, 9/10**

7/8 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid varieties
9/10 Graceys available in Standard



GRACEYS: **11/12, 15/16**

11/12 and 15/16 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid



GRACEYS: **13/14, 17/18**


13/14 Graceys available in Standard, Rigid, After Five, After Five Rigid, Mini Five and Mini Five Rigid
17/18 Graceys available in Standard and Rigid




11/14 and 12/13 Graceys are available in the Standard variety and come in a Resin 8 Colors red handle

USE OF GRACEY CURETTES


GRACEY COLOR CHART



SYSTEMATIC INSTRUMENTATION



Lingual Maxillary Left
Second Molar
SG11/12C8



Lingual Mandibular Right
First Premolar
SG5/6C8

RESIN 8 COLORS SCALERS & CURETTES



ANTERIOR SICKLE SCALERS

Purple Resin 8 Colors available with the following working ends:

Nevi 1 & Nevi 1/H5	0/00 Morse
H6/7 Hygienist	137
H5/33 Hygienist	U15/30 and U15/33 Towner - Jacquette
30/33 Jacquette	



POSTERIOR SICKLE SCALERS

Red Resin 8 Colors available with the following working ends:

204 IUFW	Nevi 2	31/32 Standard Jacquette
204S	Nevi 3	34/35 Miniature Jacquette
204SD	Nevi 4	
135		



UNIVERSAL CURETTES

Lavender Resin 8 Colors available with the following working ends:

Barnhart 1/2, 1/2R, 5/6, 5/6R	Younger-Good 7/8	Langer 1/2, 3/4, 5/6, 17/18
Columbia 2R/2L, 4R/4L, 13/14	McCall 13/14, 13/14S, 17/18	

CARING FOR YOUR RESIN INSTRUMENTS



CLEANING AND CARE OF YOUR RESIN INSTRUMENTS

In order to maintain the integrity of the resin instruments, detergents or disinfectants containing phenols or iodophors cannot be used. In addition, dry heat is not compatible with instruments with resin handles (Hu-Friedy handle #8), with resin or silicone components, inserts on any instruments, or with resin cassettes. The sterilizer equipment manufacturer's compatibility with specific materials must be observed.

EVEREDGE® TECHNOLOGY STAYS SHARPER LONGER

Revolutionary new stainless steel alloy is super-durable to stay sharper longer. That means less time sharpening and less hand fatigue.

EverEdge Technology Scalers feature a finely-honed finish to the instrument blade, for super sharp edges and efficient scaling.

PROPRIETARY PROCESS

Proprietary heat treatment and cryogenic processing ensure that the superior edge retention and wear characteristics of EVEREDGE® Technology will last the entire life of the instrument. It's not a superficial coating - EVEREDGE® Technology scalers can be sharpened again and again for your best instrument value.

RESIN 8 COLORS PART CODES

GRACEY CURETTES

GRACEY 1/2

- Standard.....SG1/2C8
- Rigid.....SG1/2RC8
- After Five®.....SRPG1/2C8
- After Five Rigid.....SRP1/2RC8
- Mini Five®.....SAS1/2C8
- Mini Five Rigid.....SS1/2RC8

GRACEY 3/4

- Standard.....SG3/4C8
- Rigid.....SG3/4RC8
- After Five.....SRPG3/4C8
- After Five Rigid.....SRP3/4RC8
- Mini Five.....SAS3/4C8
- Mini Five Rigid.....SS3/4RC8

GRACEY 5/6

- Standard.....SG5/6C8
- Rigid.....SG5/6RC8
- After Five.....SRPG5/6C8
- After Five Rigid.....SRP5/6RC8
- Mini Five.....SAS5/6C8
- Mini Five Rigid.....SS5/6RC8

GRACEY 7/8

- Standard.....SG7/8C8
- Rigid.....SG7/8RC8
- After Five.....SRPG7/8C8
- After Five Rigid.....SRP7/8RC8
- Mini Five.....SAS7/8C8
- Mini Five Rigid.....SS7/8RC8

GRACEY 9/10

- Standard.....SG9/10C8

GRACEY 11/12

- Standard.....SG11/12C8
- Rigid.....SG11/12RC8
- After Five.....SRP11/12C8
- After Five Rigid.....SP11/12RC8
- Mini Five.....SAS11/12C8
- Mini Five Rigid.....SS11/12RC8

GRACEY 13/14

- Standard.....SG13/14C8
- Rigid.....SG13/14RC8
- After Five.....SRP13/14C8
- After Five Rigid.....SP13/14RC8
- Mini Five.....SAS13/14C8
- Mini Five Rigid.....SS13/14RC8

GRACEY 15/16

- Standard.....SG15/16C8
- Rigid.....SG15/16RC8
- After Five.....SRP15/16C8
- After Five Rigid.....SP15/16RC8
- Mini Five.....SAS15/16C8
- Mini Five Rigid.....SS15/16RC8

GRACEY 17/18

- Standard.....SG17/18C8
- Rigid.....SG17/18RC8

GRACEY 11/14

- Standard.....SG11/14C8

GRACEY 12/13

- Standard.....SG12/13C8

SICKLE SCALERS

ANTERIOR

- Nevi 1.....SCNEVI1C8
- Nevi 1/H5.....SNEVI1/H5C8
- H5/33.....SH5/33C8
- H6/7.....SH6/7C8
- 30/33 Jacquette.....SJ30/33C8
- 0/00 Morse.....SM0/00C8
- 137.....SN137C8
- U15/30 Jacquette... SU15/30C8
- U15/33 Jacquette.....SU15/33C8

POSTERIOR

- 204 IUFW.....SIUFW204C8
- 204S.....S204SC8
- 204SD.....S204SDC8
- Nevi 2.....SCNEVI2C8
- Nevi 3.....SCNEVI3C8
- Nevi 4.....SCNEVI4C8
- 31/32 Jacquette.....SJ31/32C8
- 34/35 Jacquette.....SJ34/35C8
- 135.....SN135C8

UNIVERSAL CURETTES

BARNHART

- 1/2.....SBH1/2C8
- 1/2 Rigid.....SBH1/2RC8
- 5/6.....SBH5/6C8
- 5/6 Rigid.....SBH5/6RC8

COLUMBIA

- 2R/2L.....SC2R/2LC8
- 4R/4L.....SC4R/4LC8
- 13/14.....SC13/14C8

YOUNGER-GOOD

- 7/8.....SYG7/8C8

McCALL

- 13/14.....SM13/14C8
- 13S/14S.....SM13/14SC8
- 17/18.....SM17/18C8

LANGER

- 1/2.....SL1/2C8
- 3/4.....SL3/4C8
- 5/6.....SL5/6C8
- 17/18.....SL17/18C8



Denotes handle color.

Manufacturer: Hu-Friedry Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA
 Customer Care Department: Hu-Friedry Mfg. B.V. | P.O. Box 29025 | NL-3001 GA Rotterdam | Netherlands
 Free Phone 00800 48 37 43 39 | Free Fax 00800 48 37 43 40 | Website: Hu-Friedry.eu

