




## Introducing Hu-Friedy's New **ATRAUMAIR FORCEPS**

Hu-Friedy's new line of Atraumair Forceps are designed for optimal access and enhanced crown and root engagement.

Each forceps is meticulously crafted to the highest quality standard, ensuring superior performance and longevity. Atraumair Forceps are available in 13 unique patterns that provide a complete line for a variety of extraction cases.



### ATRAUMAIR POINTS OF PERFORMANCE:

- Enhanced beak geometries and longitudinal beak serrations provide a **firm, secure grasp of the crown and root reducing the risk of crushing or fracturing the tooth**
- Laser marked hinge for **greater visibility when bagged for sterilization**
- Matte finish for **enhanced contrast and reduced light reflection**
- Unique hole design reduces weight of the forceps, providing **increased comfort and control**
- Meticulously crafted from the finest surgical grade stainless steel

How the best perform

**Hu-Friedy**

# ATRAUMAIR FORCEPS

67 Forceps  
Upper Molars  
| **F67XS**



1 Standard Forceps  
Upper Incisors  
| **F1XS**



1 Apical Forceps  
Upper Incisors  
| **FAF1XS**



35 Forceps, Upper  
Canines & Premolars  
| **F35XS**



35 Apical Forceps, Upper  
Canines & Premolars  
| **FAF35XS**



51 Root Forceps  
Upper Roots, Serrated  
| **FX51XS**



36 Forceps  
Lower Incisors  
& Canines  
| **F36XS**



36 Apical Forceps  
Lower Premolars,  
Canines & Incisors  
| **FAF36XS**



13 Forceps  
Lower Premolars  
| **F13XS**



79 Forceps  
Lower Molars  
| **FX79XS**



73 Forceps  
Lower Molars  
| **FX73XS**



23 Cowhorn Forceps  
Lower Molars  
| **F23XS**



45 Root Forceps  
Lower Roots, Serrated  
| **FX45XS**