

**NEW!**

IMS SIGNATURE SERIES CASSETTES

ALWAYS PREPARED, SO YOU CAN PRACTICE MORE.

Uniquely sophisticated, ergonomic design with revolutionary new hole pattern for optimal use in Thermodisinfector.

Ergonomic latch with positive locking mechanism: safe and easy-to-use; allows for one-handed opening and visual indication if cassette is locked/unlocked.

Revolutionary new hole pattern design: with 1/3 more open area allows for improved cleaning and sterilisation. **Designed for use in thermodisinfector.**

Innovative silicone rail system features: less instrument contact for more flow through while protecting the instruments through reprocessing.

High quality, electropolished stainless steel provides a protective sheet against corrosion and a professional appearance.

3 new colors added for the silicone rail color-coding system.

New beveled edge allows for easier pick up of cassette.

Smooth Round Corners and Slotted Edges: increase drainage and reduce drying time.

WHY HU-FRIEDY IMS CASSETTES?

Hu-Friedy's Instrument Management System (IMS) is designed to keep instruments organized and intact from cleaning to chairside, providing opportunity for:

- 5-10 minutes saved per procedure
- Standardized procedures
- Increased staff safety
- Increased instrument life
- Focus on revenue-generating activities
- Reduction in overtime paid
- Potential to see 1-2 additional patients per day

How the best perform

Hu-Friedy

REFERENCE SECTION

HU-FRIEDY IMS CASSETTES

Uniquely sophisticated, ergonomic design with revolutionary new hole pattern for optimal use in TDI. IMS SIGNATURE SERIES CASSETTES are made of long-lasting, low-maintenance stainless steel. Durable, coloured silicone rails securely hold instruments throughout the processing cycle. Available in different sizes (1DIN, 1/2DIN, 1/4DIN and container cassette), configurations (Ortho, Bone Spreading, Surgical, Diagnostics and Prophylactics, etc.) and colour options, and can be specialised to meet your needs.



Cassette Color Coding

■	0 = Red	■	6 = Lavender
■	1 = Grey	■	7 = Black (new)
■	2 = White	■	8 = Blue
■	3 = Orange	■	9 = Green
■	4 = Purple (new)	■	11 = Ocean Blue (new)
■	5 = Yellow		

HU-FRIEDY IMS CASSETTES PARTCODES

Partcode	Description	Dimensions	Accessory Area	Clips
IMECC12_	Container Cassette 12 Instrument	268 x 176 x 34 mm	✓	
IMEDIN13_	Full DIN 13 Instrument	290 x 184 x 34 mm	✓	
IMEDIN18_	Full DIN 18 Instrument	290 x 184 x 34 mm	✓	
IMEDIN18_S	Full DIN 18 Instrument Surgical	290 x 184 x 34 mm	✓	✓
IMEDIN16_	Full DIN 16 Instrument	290 x 184 x 34 mm	✓	
IMEDIN20_	Full DIN 20 Instrument	290 x 184 x 34 mm		
IMEDINORT_	Full DIN Ortho (8 Instrument)	290 x 184 x 34 mm		
IMEDINOST	Full DIN Bone Spread	290 x 184 x 34 mm	✓	
IME12DIN1_	1/2 DIN 10 Instrument	152 x 184 x 34 mm		
IME12DIN8_	1/2 DIN 8 Instrument	152 x 184 x 34 mm	✓	
IME12DINOS	1/2 DIN Bone Spread	152 x 184 x 34 mm		
IME14DIN5_	1/4 DIN 5 Instrument	72 x 184 x 34 mm		

** Please refer to Hu-Friedy representative or customer care for part codes and more details.

Chicago | Milan | Rotterdam | Shanghai | Tokyo | Tuttlingen

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA
 Customer Care Department: Free Phone 00800 48 37 43 39
 Free Fax 00800 48 37 43 40 | Email: info@hu-friedy.eu
 Website: Hu-Friedy.eu



©2013 Hu-Friedy Mfg. Co., LLC. All rights reserved.

The contents of this Document including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of Hu-Friedy Mfg. Co., LLC and may not be copied, reproduced or distributed without prior written permission. All information is subject to change. IMS-GB-IME/1013