

ENZYMAX EARTH

Revision date: 15.07.2015

Revision No: 1,2

Product code: 00379-0010

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

ENZYMAX EARTH

Further trade names

IMS-1336

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Cleaning agent

Uses advised against

There is no information available on applications that are not advised.

1.3. Details of the supplier of the safety data sheet

Manufacturer

Company name: Hu-Friedy Mfg. Co., LLC

Street: 3232 N. Rockwell St.

Place: Chicago, IL 60618

Supplier

Company name: Hu-Friedy Mfg. Co., LLC.

Zweigniederlassung Deutschland

Street: Kleines Öschle 8

Place: D-78532 Tuttlingen

Telephone: 00800 48 37 43 39

Telefax: 00800 48 37 43 40

e-mail: info@hu-friedy.eu

Internet: www.hu-friedy.eu

Responsible Department: Responsible for the safety data sheet: sds@gbk-ingelheim.de

1.4. Emergency telephone number:

INTERNATIONAL: +49 - (0) 6132 - 84463, GBK GmbH (24h - 7d/w - 365d/a)

In England and Wales: NHS Direct - 0845 4647, In Scotland: NHS 24 - 08454 24 24

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture according to 1272/2008/EC

Hazard categories:

Serious eye damage/eye irritation: Eye Dam. 1

Respiratory or skin sensitisation: Resp. Sens. 1

Hazard Statements:

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Causes serious eye damage.

2.2. Label elements

Hazardous components which must be listed on the label

(Carboxymethyl)dimethyl-3-[(1-oxooctyl)amino]propylammonium hydroxide

Signal word:

Danger

Pictograms:



Hazard statements

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H318 Causes serious eye damage.

Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P260 Do not breathe vapour/spray.

ENZYMAX EARTH

Revision date: 15.07.2015

Revision No: 1,2

Product code: 00379-0010

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P310 Immediately call a POISON CENTER/doctor.
 P420 Store away from other materials.
 P501 Waste disposal according to local regulations.

2.3. Other hazards

Percentage of ingredients in admixture with unknown toxicity: 11 %

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Chemical characterization

Enzymatic cleansing agent

Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification according to Regulation (EC) No. 1272/2008 [CLP]			
73772-46-0	(Carboxymethyl)dimethyl-3-[(1-oxooctyl)amino]propylammonium hydroxide			1 - 5 %
	277-600-6			
	Eye Dam. 1; H318			
160875-66-1	2-Propylheptanoethoxilate			1 - 5 %
	Eye Dam. 1; H318			
9014-01-1	Subtilisin			0,1 - 1 %
	232-752-2	647-012-00-8	01-2119480434-38	
	Skin Irrit. 2, Eye Dam. 1, Resp. Sens. 1, STOT SE 3, Aquatic Acute 1; H315 H318 H334 H335 H400			

Full text of H and EUH phrases: see section 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

Remove contaminated soaked clothing immediately. If you feel unwell, seek medical advice.

After inhalation

Move to fresh air in case of accidental inhalation of fumes from overheating or combustion. In the event of symptoms refer for medical treatment.

After contact with skin

In case of contact with skin wash off with soap and water. Consult a doctor if skin irritation persists.

After contact with eyes

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Remove contact lens. If eye irritation persists, consult a specialist.

After ingestion

Do not provoke vomiting. Consult physician. Attention in case of vomiting - acute danger of suffocating, produced by foaming ingredients. Rinse mouth. Make drink some glasses of water. The decision whether to provoke vomiting is to be taken by a physician.

4.2. Most important symptoms and effects, both acute and delayed

Causes serious eye damage. Can cause allergic reactions.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptoms. Do not inhale vapour/aerosol. Information on label instructions.

SECTION 5: Firefighting measures

5.1. Extinguishing media

ENZYMAX EARTH

Revision date: 15.07.2015

Revision No: 1,2

Product code: 00379-0010

Suitable extinguishing media

Product does not burn, fire-extinguishing activities according to surrounding.

Unsuitable extinguishing media

Full water jet.

5.2. Special hazards arising from the substance or mixture

Fire may produce: Carbon monoxide (CO), carbon dioxide (CO₂) and nitrogen oxides (NO_x).

5.3. Advice for firefighters

Use breathing apparatus with independent air supply.

Protective suit.

Additional information

Cool containers at risk with water spray jet. Fire residues and contaminated firefighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Keep away noninvolved persons. Avoid contact with skin, eyes and clothing. Ensure adequate ventilation.

6.2. Environmental precautions

Do not discharge into the drains/surface waters/ground water.

6.3. Methods and material for containment and cleaning up

Soak up with inert absorbent material (e.g. sand, silica gel, acid binder, universal binder).

6.4. Reference to other sections

Information for disposal see section 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid contact with eyes, skin or mucous membrane. Do not inhale vapour/aerosol.

Advice on protection against fire and explosion

No special protective measures against fire required.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed in a dry, cool and well-ventilated place.

Advice on storage compatibility

Incompatible with strong acids and oxidizing agents.

Further information on storage conditions

Keep away from food, drink and animal feeding stuffs.

7.3. Specific end use(s)

No data available.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Exposure limits (EH40)

CAS No	Substance	ppm	mg/m ³	fibres/ml	Category	Origin
57-55-6	Propane-1,2-diol, particulates	-	10		TWA (8 h)	WEL
		-	-		STEL (15 min)	WEL
9014-01-1	Subtilisins (Bacillus subtilis Carlsberg)	-	0.00004		TWA (8 h)	WEL
		-	-		STEL (15 min)	WEL

8.2. Exposure controls

ENZYMAX EARTH

Revision date: 15.07.2015

Revision No: 1,2

Product code: 00379-0010

Appropriate engineering controls

Ensure adequate ventilation, especially in confined areas.

Protective and hygiene measures

Wash hands before breaks and immediately after handling the product. When using do not eat, drink or smoke. Take off immediately all contaminated clothing. Avoid contact with skin, eyes and clothing.

Eye/face protection

Tightly fitting goggles (EN 166).

Eye wash bottle with pure water (EN 15154).

Hand protection

Protective gloves (EN 374); Rubber gloves. (recommended: minimum protection index 2, corresponding to a permeation rate > 30 minutes according to EN 374)

Requirements can vary as a function of the use. Therefore it is necessary to adhere additionally to the recommendations given by the manufacturer of protective gloves.

Skin protection

Long sleeved clothing (EN 368).

Respiratory protection

No personal respiratory protective equipment normally required.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state:	Liquid
Colour:	Green
Odour:	characteristic

pH-Value:	6,5 - 7,5
-----------	-----------

Changes in the physical state

Flash point:	n.a.
Lower explosion limits:	n.d.
Upper explosion limits:	
Ignition temperature:	n.a.
Density:	1,10 - 1,11 g/cm ³
Water solubility: (at 20 °C)	Completely miscible

9.2. Other information

No data available.

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Keep away from heat.

10.4. Conditions to avoid

To avoid thermal decomposition, do not overheat.

10.5. Incompatible materials

Strong oxidizing agents.

10.6. Hazardous decomposition products

Carbon monoxide (CO), carbon dioxide (CO₂) and nitrogen oxides (NO_x).

Further information

No decomposition if stored and applied as directed.

ENZYMAX EARTH

Revision date: 15.07.2015

Revision No: 1,2

Product code: 00379-0010

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

Based on available data, the classification criteria are not met.

Irritation and corrosivity

Causes serious eye damage.

Sensitising effects

May cause allergy or asthma symptoms or breathing difficulties if inhaled. (Subtilisin)

STOT-single exposure

Based on available data, the classification criteria are not met.

Severe effects after repeated or prolonged exposure

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Practical experience

Other observations

Contact with eyes, skin or oral tissues may cause irritation. Ingestion causes irritation of upper respiratory system and gastrointestinal disturbance. Repeated or prolonged skin contact may cause skin irritation and/or dermatitis and sensitization of susceptible persons. May cause sensitization of susceptible persons by inhalation.

SECTION 12: Ecological information

12.1. Toxicity

No data available.

12.2. Persistence and degradability

The surfactant(s) contained in this preparation complies(comply) with the biodegradability criteria as laid down in Regulation (EC) No.648/2004 on detergents. Data to support this assertion are held at the disposal of the competent authorities of the Member States and will be made available to them, at their direct request or at the request of a detergent manufacturer.

12.3. Bioaccumulative potential

No data available.

12.4. Mobility in soil

No data available.

12.5. Results of PBT and vPvB assessment

No data available.

12.6. Other adverse effects

No data available.

Further information

Do not flush into surface water or sanitary sewer system.

Low hazard to waters.

Ecological data are not available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Advice on disposal

Can be disposed of as a solid waste or burned in a suitable installation subject to local regulations.

Where possible recycling is preferred to disposal.

Waste disposal number of waste from residues/unused products

070699

WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of fats, grease, soaps, detergents, disinfectants and cosmetics; wastes not otherwise specified

ENZYMAX EARTH

Revision date: 15.07.2015

Revision No: 1,2

Product code: 00379-0010

Waste disposal number of contaminated packaging

150102 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); plastic packaging

Contaminated packaging

Empty containers should be taken for local recycling, recovery or waste disposal.

Contaminated packaging should be emptied as far as possible and after appropriate cleansing may be taken for reuse.

Packaging that cannot be cleaned should be disposed of like the product.

SECTION 14: Transport information

Land transport (ADR/RID); Marine transport (IMDG); Air transport (ICAO); Inland waterways transport (ADN):

14.1. UN number:

Not applicable. No hazardous material as defined by the transport regulations.

14.2. UN proper shipping name:

Not applicable. No hazardous material as defined by the transport regulations.

14.3. Transport hazard class(es):

Not applicable. No hazardous material as defined by the transport regulations.

14.4. Packing group:

Not applicable. No hazardous material as defined by the transport regulations.

14.5. Environmental hazards

Not applicable. No hazardous material as defined by the transport regulations.

14.6. Special precautions for user

Not applicable. No hazardous material as defined by the transport regulations.

14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Not applicable. No hazardous material as defined by the transport regulations.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU regulatory information

2004/42/EC (VOC): 1 %

National regulatory information

Employment restrictions: Observe employment restrictions for young people. Observe employment restrictions for child bearing mothers and nursing.

Water contaminating class (D): 1 - slightly water contaminating

15.2. Chemical safety assessment

For this substance a chemical safety assessment is not required.

SECTION 16: Other information

Changes

Changes in chapter: 1 - 16

Abbreviations and acronyms

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route

RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses

ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure

IMDG = International Maritime Code for Dangerous Goods

IATA/ICAO = International Air Transport Association / International Civil Aviation Organization

MARPOL = International Convention for the Prevention of Pollution from Ships

IBC-Code = International Code for the Construction and Equipment of Ships Carrying Dangerous Chemicals in Bulk

ENZYMAX EARTH

Revision date: 15.07.2015

Revision No: 1,2

Product code: 00379-0010

GHS = Globally Harmonized System of Classification and Labelling of Chemicals

REACH = Registration, Evaluation, Authorization and Restriction of Chemicals

CAS = Chemical Abstract Service

EN = European norm

ISO = International Organization for Standardization

DIN = Deutsche Industrie Norm

PBT = Persistent Bioaccumulative and Toxic

vPvB = Very Persistent and very Bio-accumulative

LD = Lethal dose

LC = Lethal concentration

EC = Effect concentration

IC = Median immobilisation concentration or median inhibitory concentration

Relevant H- and EUH-phrases (Number and full text)

H315 Causes skin irritation.

H318 Causes serious eye damage.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled .

H335 May cause respiratory irritation.

H400 Very toxic to aquatic life.

Further Information

Data of items 4 to 8, as well as 10 to 12, do partly not refer to the use and the regular employing of the product (in this sense consult information on use and on product), but to liberation of major amounts in case of accidents and irregularities.

The information describes exclusively the safety requirements for the product(s) and is based on the present level of our knowledge.

The delivery specifications are contained in the corresponding product sheet.

This data does not constitute a guarantee for the characteristics of the product(s) as defined by the legal warranty regulations.

(n.a. = not applicable; n.d. = not determined)

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)